

Tbx15 Antibody, FITC conjugated



PACO30296

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable transcriptional regulator involved in the development of the skeleton of the limb, vertebral column and head. Acts by controlling the number of mesenchymal precursor cells and chondrocytes.

Gene ID:

Tbx15

Uniprot

O70306

Synonyms:

T-box transcription factor TBX15 (T-box protein 15) (MmTBx8) (T-box transcription factor TBX14) (T-box protein 14), Tbx15, Tbx14 Tbx8

Immunogen:

Recombinant Mouse T-box transcription factor TBX15 protein (1-120AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A