

# Tbx15 Antibody, Biotin conjugated



PACO30297

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probable transcriptional regulator involved in the development of the skeleton of the limb, vertebral column and head. Acts by controlling the number of mesenchymal precursor cells and chondrocytes.

**Gene ID:**

Tbx15

**Uniprot**

O70306

**Synonyms:**

T-box transcription factor TBX15 (T-box protein 15) (MmTBx8) (T-box transcription factor TBX14) (T-box protein 14), Tbx15, Tbx14 Tbx8

**Immunogen:**

Recombinant Mouse T-box transcription factor TBX15 protein (1-120AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A