

# TCN1 Antibody, Biotin conjugated



PACO30313

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds vitamin B12 with femtomolar affinity and protects it from the acid, c environment of the stomach.

**Gene ID:**

TCN1

**Uniprot**

P20061

**Synonyms:**

Transcobalamin-1 (TC-1) (Haptocorrin) (HC) (Protein R) (Transcobalamin I) (TC I) (TCI), TCN1, TC1

**Immunogen:**

Recombinant Human Transcobalamin-1 protein (24-433AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A