

# TCTA Antibody, HRP conjugated



PACO30319

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be required for cellular fusion during osteoclastogenesis.

**Gene ID:**

TCTA

**Uniprot**

P57738

**Synonyms:**

T-cell leukemia translocation-altered gene protein (T-cell leukemia translocation-associated gene protein), TCTA

**Immunogen:**

Recombinant Human T-cell leukemia translocation-altered gene protein (62-103AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A