

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This protein has been shown to exist at transcript level. Moderate nuclear positivity was observed in cells in seminiferous ducts of testis. Moderate cytoplasmic positivity was observed in a subset of lymphoid cells outside reaction centers, glial cells, respiratory epithelia and ciliated cells in fallopian tube. Neuronal cells and cells in molecular layer of cerebellum were strongly stained. Remaining normal tissues were in general negative.

Gene ID:

TEX12

Uniprot

Q9BXU0

Synonyms:

Testis-expressed protein 12, TEX12

Immunogen:

Recombinant Human Testis-expressed protein 12 protein (1-123AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A