

Thbs2 Antibody



PACO30374

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Thbs2

Uniprot

Q03350

Synonyms:

Thrombospondin-2, Thbs2, Tsp2

Immunogen:

Recombinant Mouse Thrombospondin-2 protein (19-232AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A