

TIMM10 Antibody, FITC conjugated



PACO30400

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Mitochondrial intermembrane chaperone that participates in the import and insertion of multi-pass transmembrane proteins into the mitochondrial inner membrane. May also be required for the transfer of beta-barrel precursors from the TOM complex to the sorting and assembly machinery (SAM complex) of the outer membrane. Acts as a chaperone-like protein that protects the hydrophobic precursors from aggregation and guide them through the mitochondrial intermembrane space.

Gene ID:

TIMM10

Uniprot

P62072

Synonyms:

Mitochondrial import inner membrane translocase subunit Tim10, TIMM10, TIM10

Immunogen:

Recombinant Human Mitochondrial import inner membrane translocase subunit Tim10 protein (1-90AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A