

PACO30410

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be a allosteric regulator of serine/threonine-protein phosphatase 2A (PP2A). Isoform 1 inhibits catalytic activity of the PP2A(D) core complex in vitro. The PP2A(C):TIPRL complex does not show phosphatase activity. May play a role in the regulation of ATM/ATR signaling pathway controlling DNA replication and repair.

**Gene ID:**

TIPRL

**Uniprot**

O75663

**Synonyms:**

TIP41-like protein (Putative MAPK-activating protein PM10) (Type 2A-interacting protein) (TIP), TIPRL

**Immunogen:**

Recombinant Human TIP41-like protein (1-178AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A