

TM4SF5 Antibody



PACO30438

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

integral component of plasma membrane.

Gene ID:

TM4SF5

Uniprot

O14894

Synonyms:

Transmembrane 4 L6 family member 5 (Tetraspan transmembrane protein L6H), TM4SF5


Immunogen:

Recombinant Human Transmembrane 4 L6 family member 5 protein (112-157AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO30438 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO30438 at dilution of 1:100.