

TM7SF3 Antibody



PACO30442

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

extracellular exosome.

Gene ID:

TM7SF3

Uniprot

Q9NS93

Synonyms:

Transmembrane 7 superfamily member 3 (Seven span transmembrane protein),
TM7SF3


Immunogen:

Recombinant Human Transmembrane 7 superfamily member 3 protein (22-295AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO30442 at dilution of 1:100.



Immunofluorescent analysis of U251 cells using PACO30442 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).