

TMEM115 Antibody



PACO30462

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May play a role in retrograde transport of proteins from the Golgi to the endoplasmic reticulum. May indirectly play a role in protein glycosylation in the Golgi.

Gene ID:

TMEM115

Uniprot

Q12893

Synonyms:

Transmembrane protein 115 (Placental protein 6) (Protein PL6), TMEM115, PL6


Immunogen:

Recombinant Human Transmembrane protein 115 protein (187-351AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO30462 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO30462 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human testis tissue using PACO30462 at dilution of 1:100.