

TMEM177 Antibody



PACO30518

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200, IF:1:50-1:200

Protein Background:

Gene ID:

TMEM177

Uniprot

Q53S58

Synonyms:

Transmembrane protein 177, TMEM177


Immunogen:

Recombinant Human Transmembrane protein 177 protein (219-311AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO30518 at dilution of 1:100.



Immunofluorescent analysis of PC-3 cells using PACO30518 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human colon cancer using PACO30518 at dilution of 1:100.