

TMEM184C Antibody



PACO30530

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Possible tumor suppressor which may play a role in cell growth.

Gene ID:

TMEM184C

Uniprot

Q9NVA4

Synonyms:

Transmembrane protein 184C (Transmembrane protein 34), TMEM184C, TMEM34

Immunogen:

Recombinant Human Transmembrane protein 184C protein (107-178AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO30530 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO30530 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human thymus tissue using PACO30530 at dilution of 1:100.