

TMEM186 Antibody



PACO30534

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200, IF:1:50-1:200

Protein Background:

mitochondrion.

Gene ID:

TMEM186

Uniprot

Q96B77

Synonyms:

Transmembrane protein 186, TMEM186, C16orf51

Immunogen:

Recombinant Human Transmembrane protein 186 protein (124-213AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO30534 at dilution of 1:100.



Immunofluorescent analysis of MCF-7 cells using PACO30534 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO30534 at dilution of 1:100.