

TMEM196 Antibody



PACO30546

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TMEM196

Uniprot

Q5HYL7

Synonyms:

Transmembrane protein 196, TMEM196

Immunogen:

Recombinant Human Transmembrane protein 196 protein (53-89AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A