

TMEM205 Antibody



PACO30550

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

In cancer cells, plays a role in resistance to the chemotherapeutic agent cisplatin.

Gene ID:

TMEM205

Uniprot

Q6UW68

Synonyms:

Transmembrane protein 205, TMEM205


Immunogen:

Recombinant Human Transmembrane protein 205 protein (104-153AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO30550 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO30550 at dilution of 1:100.