

TMEM217 Antibody




PACO30570

Product Information

| | |
|--|--|
| Size: | Protein Background: |
| 50ug | |
| Reactivity: | Gene ID: |
| Human | TMEM217 |
| Source: | Uniprot |
| Rabbit | Q8N7C4 |
| Isotype: | Synonyms: |
| IgG | Transmembrane protein 217, TMEM217, C6orf128 |
| Applications: | Immunogen: |
| ELISA, WB, IF | Recombinant Human Transmembrane protein 217 protein (150-211AA). |
| Recommended dilutions: | Storage: |
| ELISA:1:2000-1:10000, WB:1:1000-1:5000, IF:1:50-1:200 | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |

Product Images



Western blot

All lanes: TMEM217 antibody at 10 μ g/ml + Jurkat whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 27, 25 kDa

Observed band size: 27 kDa

.



Immunofluorescent analysis of HepG2 cells using PACO30570 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).