

# TMEM98 Antibody



PACO30670

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TMEM98

**Uniprot**

Q9Y2Y6

**Synonyms:**

Transmembrane protein 98 (Protein TADA1), TMEM98

**Immunogen:**

Recombinant Human Transmembrane protein 98 protein (25-226AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A