

Tnnt2 Antibody, FITC conjugated



PACO30720

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Troponin T is the tropomyosin-binding subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

Gene ID:

Tnnt2

Uniprot

P50752

Synonyms:

Troponin T, cardiac muscle (TnTc) (Cardiac muscle troponin T) (cTnT), Tnnt2

Immunogen:

Recombinant Mouse Troponin T, cardiac muscle protein (2-301AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A