

TPPP2 Antibody



PACO30734

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200, IF:1:50-1:200

Protein Background:

May bind tubulin but has no microtubule bundling activity.

Gene ID:

TPPP2

Uniprot

P59282

Synonyms:

Tubulin polymerization-promoting protein family member 2 (Protein p25-beta)
(TPPP/p18), TPPP2, C14orf8


Immunogen:

Recombinant Human Tubulin polymerization-promoting protein family member 2
protein (1-170AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: TPPP2 antibody at 8 μ g/ml + HepG2 whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution


Predicted band size: 19 kDa

Observed band size: 19 kDa

.



Immunofluorescent analysis of MCF-7 cells using PACO30734 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO30734 at dilution of 1:100.