

Tpsb2 Antibody, Biotin conjugated



PACO30757

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Tryptase is the major neutral protease present in mast cells and is secreted upon the coupled activation-degranulation response of this cell type. Has an immunoprotective role during bacterial infection. Required to efficiently combat K. pneumoniae infection.

Gene ID:

Tpsb2

Uniprot

P50343

Synonyms:

Tryptase beta-2 (Tryptase-2) (EC 3.4.21.59) (Mast cell protease 6) (rMCP-6), Tpsb2, Mcp6
Mcp6

Immunogen:

Recombinant Rat Tryptase beta-2 protein (30-274AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A