

PACO30758

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Catalyzes the O-sulfation of tyrosine residues within acid, c motifs of polypeptides.

**Gene ID:**

TPST2

**Uniprot**

O60704

**Synonyms:**

Protein-tyrosine sulfotransferase 2 (EC 2.8.2.20) (Tyrosylprotein sulfotransferase 2) (TPST-2), TPST2

**Immunogen:**


Recombinant Human Protein-tyrosine sulfotransferase 2 protein (26-377AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO30758 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO30758 at dilution of 1:100.