

Tpt1 Antibody



PACO30762

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:10-1:100

Protein Background:

Involved in calcium binding and microtubule stabilization.

Gene ID:

Tpt1

Uniprot

P63028

Synonyms:

Translationally-controlled tumor protein (TCTP) (21 kDa polypeptide) (p21) (p23), Tpt1, Trt


Immunogen:

Recombinant Mouse Translationally-controlled tumor protein (1-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO30762 at dilution of 1:20.