

TPTE Antibody, HRP conjugated



PACO30767

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Could be involved in signal transduction.

Gene ID:

TPTE

Uniprot

P56180

Synonyms:

Putative tyrosine-protein phosphatase TPTE (EC 3.1.3.48) (Cancer/testis antigen 44) (CT44) (Transmembrane phosphatase with tensin homology) (Tumor antigen BJ-HCC-5), TPTE

Immunogen:

Recombinant Human Putative tyrosine-protein phosphatase TPTE protein (175-380AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A