

TRIM16 Antibody



PACO30802

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in the regulation of keratinocyte differentiation.

Gene ID:

TRIM16

Uniprot

O95361

Synonyms:

Tripartite motif-containing protein 16 (Estrogen-responsive B box protein), TRIM16, EBBP

Immunogen:

Recombinant Human Tripartite motif-containing protein 16 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A