

TSPAN8 Antibody, FITC conjugated



PACO30864

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

extracellular exosome, integral component of plasma membrane, cell surface receptor signaling pathway.

Gene ID:

TSPAN8

Uniprot

P19075

Synonyms:

Tetraspanin-8 (Tspan-8) (Transmembrane 4 superfamily member 3) (Tumor-associated antigen CO-029), TSPAN8, TM4SF3

Immunogen:

Recombinant Human Tetraspanin-8 protein (110-205AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A