

TTL Antibody, FITC conjugated



PACO30868

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the post-translational addition of a tyrosine to the C-terminal end of dephosphorylated alpha-tubulin.

Gene ID:

TTL

Uniprot

Q8NG68

Synonyms:

Tubulin--tyrosine ligase (TTL) (EC 6.3.2.25), TTL

Immunogen:

Recombinant Human Tubulin--tyrosine ligase protein (1-377AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A