

Ttr Antibody, FITC conjugated



PACO30882

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Gene ID:

Ttr

Uniprot

P07309

Synonyms:

Transthyretin (Prealbumin), Ttr

Immunogen:

Recombinant Mouse Transthyretin protein (23-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A