

tuf1 Antibody



PACO30896

Product Information

Size:

50ug

Reactivity:

Acinetobacter baumannii

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This protein promotes the GTP-dependent binding of aminoacyl-tRNA to the A-site of ribosomes during protein biosynthesis.

Gene ID:

tuf1

Uniprot

A3M1F6

Synonyms:

Elongation factor Tu (EF-Tu), tuf1; tuf2

Immunogen:

Recombinant Acinetobacter baumannii Elongation factor Tu protein (1-396AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A