

PACO30904

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Possesses significant protein thiol-disulfide oxidase activity.

**Gene ID:**

TXNDC12

**Uniprot**

O95881

**Synonyms:**

Thioredoxin domain-containing protein 12 (EC 1.8.4.2) (Endoplasmic reticulum resident protein 18) (ER protein 18) (ERp18) (Endoplasmic reticulum resident protein 19) (ER protein 19) (ERp19) (Thioredoxin-like protein p19) (hTLP19), TXNDC12, TLP19

**Immunogen:**

Recombinant Human Thioredoxin domain-containing protein 12 protein (27-172AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A