

PACO30968

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in endoplasmic reticulum-associated degradation (ERAD) for misfolded luminal proteins, possibly by tethering VCP to the endoplasmic reticulum membrane. May play a role in reproduction.

Gene ID:

UBXN8

Uniprot

O00124

Synonyms:

UBX domain-containing protein 8 (Reproduction 8 protein) (Rep-8 protein) (UBX domain-containing protein 6), UBXN8, D8S2298E REP8 UBXD6


Immunogen:

Recombinant Human UBX domain-containing protein 8 protein (55-270AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO30968 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO30968 at dilution of 1:100.