

RPS10 Antibody, Biotin conjugated



PACO31003

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of the 40S ribosomal subunit.

Gene ID:

RPS10

Uniprot

P46783

Synonyms:

40S ribosomal protein S10 (Small ribosomal subunit protein eS10), RPS10

Immunogen:

Recombinant Human 40S ribosomal protein S10 protein (4-165AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A