

VAMP1 Antibody, HRP conjugated



PACO31029

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the targeting and/or fusion of transport vesicles to their target membrane.

Gene ID:

VAMP1

Uniprot

P23763

Synonyms:

Vesicle-associated membrane protein 1 (VAMP-1) (Synaptobrevin-1), VAMP1, SYB1

Immunogen:

Recombinant Human Vesicle-associated membrane protein 1 protein (1-96AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A