

# VAMP1 Antibody, FITC conjugated



PACO31030

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in the targeting and/or fusion of transport vesicles to their target membrane.

**Gene ID:**

VAMP1

**Uniprot**

P23763

**Synonyms:**

Vesicle-associated membrane protein 1 (VAMP-1) (Synaptobrevin-1), VAMP1, SYB1

**Immunogen:**

Recombinant Human Vesicle-associated membrane protein 1 protein (1-96AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A