

# WNT4 Antibody, HRP conjugated



PACO31093

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.

Overexpression may be associated with abnormal proliferation in human breast tissue.

**Gene ID:**

WNT4

**Uniprot**

P56705

**Synonyms:**

Protein Wnt-4, WNT4

**Immunogen:**

Recombinant Human Protein Wnt-4 protein (23-351AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A