

YPEL3 Antibody, Biotin conjugated



PACO31131

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in proliferation and apoptosis in myeloid precursor cells.

Gene ID:

YPEL3

Uniprot

P61236

Synonyms:

Protein yippee-like 3, YPEL3

Immunogen:

Recombinant Human Protein yippee-like 3 protein (1-119AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A