

ZNF131 Antibody



PACO31152

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Plays a role during development and organogenesis as well as in the function of the adult central nervous system. May be involved in transcriptional regulation as a repressor of ESR1/ER-alpha signaling.

Gene ID:

ZNF131

Uniprot

P52739

Synonyms:

Zinc finger protein 131, ZNF131


Immunogen:

Recombinant Human Zinc finger protein 131 protein (80-379AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31152 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human thymus tissue using PACO31152 at dilution of 1:100.