

ZNF131 Antibody, HRP conjugated



PACO31153

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role during development and organogenesis as well as in the function of the adult central nervous system. May be involved in transcriptional regulation as a repressor of ESR1/ER-alpha signaling.

Gene ID:

ZNF131

Uniprot

P52739

Synonyms:

Zinc finger protein 131, ZNF131

Immunogen:

Recombinant Human Zinc finger protein 131 protein (80-379AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A