

PACO31168

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZKSCAN7

**Uniprot**

Q9POL1

**Synonyms:**

Zinc finger protein with KRAB and SCAN domains 7 (Zinc finger protein 167) (Zinc finger protein 448) (Zinc finger protein 64), ZKSCAN7, ZNF167 ZNF448 ZNF64

**Immunogen:**

Recombinant Human Zinc finger protein with KRAB and SCAN domains 7 protein (1-276AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A