

# ZNF175 Antibody



PACO31172

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Down-regulates the expression of several chemokine receptors. Interferes with HIV-1 replication by suppressing Tat-induced viral LTR promoter activity.

**Gene ID:**

ZNF175

**Uniprot**

Q9Y473

**Synonyms:**

Zinc finger protein 175 (Zinc finger protein OTK18), ZNF175

**Immunogen:**


Recombinant Human Zinc finger protein 175 protein (1-711AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO31172 at dilution of 1:100.