

ZNF230 Antibody



PACO31224

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF230

Uniprot

Q9UIE0

Synonyms:

Zinc finger protein 230 (Zinc finger protein FDZF2), ZNF230, FDZF2


Immunogen:

Recombinant Human Zinc finger protein 230 protein (201-474AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO31224 at dilution of 1:100.