

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcription factor that acts as a cardiac regulator and an effector of alpha1-adrenergic signaling. Binds to PE response elements (PERE) present in the promoter of genes such as ANF/NPPA and acts as a direct transcriptional activator of NPPA. Also acts as a cofactor with GATA4, a key cardiac regulator.

Gene ID:

ZNF260

Uniprot

Q3ZCT1

Synonyms:

Zinc finger protein 260 (Zfp-260), ZNF260, ZFP260


Immunogen:

Recombinant Human Zinc finger protein 260 protein (1-200AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skin tissue using PACO31240 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO31240 at dilution of 1:100.