

ZNF300 Antibody



PACO31256

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Has a transcriptional repressor activity.

Gene ID:

ZNF300

Uniprot

Q96RE9

Synonyms:

Zinc finger protein 300, ZNF300

Immunogen:

Recombinant Human Zinc finger protein 300 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A