

ZNF440 Antibody



PACO31316

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF440

Uniprot

Q8IYI8

Synonyms:

Zinc finger protein 440, ZNF440

Immunogen:

Recombinant Human Zinc finger protein 440 protein (1-280AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A