

# ZNF484 Antibody



PACO31328

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF484

**Uniprot**

Q5JVG2

**Synonyms:**

Zinc finger protein 484, ZNF484

**Immunogen:**


Recombinant Human Zinc finger protein 484 protein (118-318AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO31328 at dilution of 1:100.