

# ZNF486 Antibody



PACO31332

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZNF486

**Uniprot**

Q96H40

**Synonyms:**

Zinc finger protein 486 (KRAB domain only protein 2), ZNF486, KRBO2

**Immunogen:**

Recombinant Human Zinc finger protein 486 protein (251-463AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A