

ZNF546 Antibody



PACO31360

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF546

Uniprot

Q86UE3

Synonyms:

Zinc finger protein 546 (Zinc finger protein 49), ZNF546, ZNF49


Immunogen:

Recombinant Human Zinc finger protein 546 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: ZNF546 antibody at 4 μ g/ml

Lane 1: A549 whole cell lysate

Lane 2: 293T whole cell lysate

Lane 3: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 99 kDa

Observed band size: 99 kDa