

ZNF563 Antibody, FITC conjugated



PACO31386

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in transcriptional regulation.

Gene ID:

ZNF563

Uniprot

Q8TA94

Synonyms:

Zinc finger protein 563, ZNF563

Immunogen:

Recombinant Human Zinc finger protein 563 protein (214-476AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A