

ZNF582 Antibody, HRP conjugated



PACO31393

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in transcriptional regulation.

Gene ID:

ZNF582

Uniprot

Q96NG8

Synonyms:

Zinc finger protein 582, ZNF582

Immunogen:

Recombinant Human Zinc finger protein 582 protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A