

ZNF583 Antibody



PACO31396

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF583

Uniprot

Q96ND8

Synonyms:

Zinc finger protein 583 (Zinc finger protein L3-5), ZNF583


Immunogen:

Recombinant Human Zinc finger protein 583 protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO31396 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO31396 at dilution of 1:100.